

pE-4000

UNIVERSAL LIGHT SOURCE



ALL WAVELENGTHS INCLUDED

• POWERFUL • EFFICIENT • COMPACT



www.CoolLED.com

CoolLED 

CoolLED pE-4000 – The Universal Light Source

The pE-4000 sets the standard as the universal light source for fluorescence microscopy. Users can operate the system as a simple white light source (replacing a conventional mercury-based light source), or as an advanced, fully-controllable, excitation and stimulation source. The flexibility and extensive functionality of the pE-4000 broadens the range of illumination options in core facilities.

At the centre of the pE-4000 is CoolLED's novel, patent pending, wavelength-grouping concept which offers more power in an efficient system design. Wavelength-grouping ensures optimal compatibility with all single and multi-band filter sets.

WHITE for Simplicity



WHITE

ADVANCED for Control



ADVANCED




IDEAL FOR MULTI-USER AND CORE FACILITIES

- Simple on/off
- Precise intensity control
- Easy to use, no training required
- Pre-sets allowing lab manager to match white spectrum to existing filter cubes
- Higher contrast images from matched-white spectrum


pE-4000 LED SYSTEM

- Excellent field uniformity at sample
- No mercury
- Long Life: 25,000 hours
- No bulb changing, bulb alignment or warm up process
- Quiet operation
- High efficiency
- Wide range of microscope adaptors




IDEAL FOR ADVANCED RESEARCH

- Individual LED wavelength selection
- Rapid switching between LED wavelengths enables capture of high speed events
- TTL & USB interfaces with imaging packages
- Excitation filters can be fitted in optical path for controlled switching with no moving parts
- Analogue for dynamic intensity control
- Optical feedback for applications requiring higher stability
- Internal function generator for electrophysiology and optogenetics applications
- Compatible with all single and multi-band filter sets



16 SELECTABLE WAVELENGTHS NO MODULARITY • COMPATIBLE WITH ALL FILTER SETS



BROADEST SPECTRUM • BRIGHTEST LEDs



NO MORE MERCURY!



SPECIFICATIONS

- pE-4000 system:** Main unit with complete set of wavelengths, dual function manual control pod, power supply block, mains cable, USB cable
- Light Delivery:** Single liquid light guide or fiber options
- Collimating optics:** Range of collimating optics and adaptors to fit most fluorescence microscopes
- LED Wavelengths:** Wavelengths are divided across 4 channels with each channel having individual control. Power measured at sample plane.

CHANNEL 1			CHANNEL 2			CHANNEL 3			CHANNEL 4		
Wavelength (nm)	Power (mW)	FWHM (nm)	Wavelength (nm)	Power (mW)	FWHM (nm)	Wavelength (nm)	Power (mW)	FWHM (nm)	Wavelength (nm)	Power (mW)	FWHM (nm)
365	15.12	12	460	240.19	19	525	39.00	34	635	123.88	16
385	59.39	11	470	192.29	20	550	166.74	85	660	159.01	21
405	188.32	18	490	56.00	32	580	166.74	85	740	36.62	29
435	138.75	16	500	59.26	29	595	32.52	15	770	15.61	28

Powers measured at sample plane of research grade microscope using 100% mirror in cube and 10X objectives. (Figures are typical values). Note wavelength settings 550nm and 580nm use single broad, high intensity peak (see spectrum on pages 2-3)

CONTROL & INTERFACE

- Manual:** Dual function manual control pod for White mode or Advanced mode
- Remote:** Via USB for independent on/off and intensity control of each channel. Triggering speed <1ms
Via 4 TTL inputs for independent on/off control of each channel. Triggering speed <20us
Via single TTL for on/off control of manual or software selected channels
Via 4 analogue inputs 0-10V, 0-300kHz for dynamic control of intensity from external analogue signals
- Synch Out:** 4 TTL outputs for each channel – active high
1 TTL output for any channel – active high
- Programmable interface:** 4 TTL outputs for on/off control of peripherals (transmitted light sources, stages etc)
4 analogue outputs for intensity control of peripherals (can be programmed to mirror LED intensities for channel control) 0-10V full scale.
- Function Generator:** Internally generated sine, pulse and ramps for each channel programmed via pod.
- Connectivity:** USB (B type) for PC connection. All other TTL and Analogue inputs/outputs via 25way 'D-type' female connector (optional rear mounting expansion box available for BNC connectivity).
- Imaging Software:** Recognised as 'CoolLED pE-2 peripheral' under common software e.g. Micromanager, MetaMorph, cellSens, NIS Elements, ImagePro, etc.

POWER

- Power requirements:** 110-240Va.c. 50/60Hz, 2.5A
- Power consumption:**
- | | |
|-----------------------------|----------|
| Standby (i.e. no LEDs on) | Max 7W |
| Single wavelength operation | Max 41W |
| Dual wavelength operation | Max 75W |
| Triple wavelength operation | Max 93W |
| Quad wavelength operation | Max 112W |

DIMENSIONS

- Main unit:** 200mm(w) x 75mm(d) x 185mm(h) – Weight 3.5kg
- Control pod:** 154mm(w) x 135mm(d) x 40mm(h) – Weight 0.95kg
- Power Supply:** 164mm(w) x 64mm(d) x 35mm(h) – Weight 0.58kg

TO ORDER

- pE-4000-L-SYS-ZZ** Main unit, control pod, power supply plus cables for use with 3mm liquid light guide
- pE-4000-F-SYS-ZZ** Main unit, control pod, power supply plus cables for use with SMA terminated fiber
- pE-1904** 3mm diameter, 1m long liquid light guide
- pE-1908** 3mm diameter, 3m long liquid light guide
- pE-10400-YYY** Microscope adaptor with collimating optics. To specify microscope code (YYY) see <http://www.cooled.com/Life-Sciences-Analytical/Products/Microscope-Adaptors/> for further information
- Fiber Options** see <http://www.cooled.com/Life-Sciences-Analytical/Products/Accessories/Light-Delivery/> for further information
- pE-4000-EB25D** Rear mounting expansion box for 25way D-type to BNC connectivity
- Specify local power cable (ZZ). 10=Australia, 20=Europe, 30=UK, 40=USA
- Warranty** System: 12 months, LEDs: 25,000 hours

ENVIRONMENT & SAFETY

LED products are more sustainable and energy efficient than conventional light sources. CoolLED's products have the following benefits:

- No Mercury
- Energy Efficient: 80% less power
- Long lifetime
- No bulb replacements
- Reduced risk of eye damage
- No special disposal regulations



CONTACT



- Online:** www.cooled.com
- Phone:** +44 (0) 1264 323040 (Worldwide)
1-800-877-0128 (USA + Canada)
- Email:** info@cooled.com

