

JuniorLED

HIGHLIGHTS

- Entry-level microscopes for brightfield observation
- Integrated X-Y mechanical stage
- Double LED illumination
- Cordless operation
- Ergonomic carrying grip
- 2 Years warranty

TECHNICAL SPECIFICATIONS

EYEPIECE

- Wide field WF 10x/16 mm eyepiece
- Secured and with pointer

HEAD

- Monocular models with 45° inclined tubes
- 360° rotatable and secured

NOSEPIECE

Revolving nosepiece for 3 objectives on ball bearings

OBJECTIVES

- Achromatic 4x/0.10, 10x/0.25 and S40x/0.65 or Achromatic 4x/0.10, 10x/0.25 and S60x/0.85
 - The S40x, S60x and S100x oil immersion objectives are spring loaded
- All optics are anti-fungus treated and anti-reflection coated for maximum light throughput*

FOCUSING

- Coaxial coarse and fine adjustments, 90 graduations, 0.02 µm per graduation, 2.35 mm per rotation, total travel approximately 20 mm
- Supplied with an adjustable rack stop to prevent damage to sample and objectives

STAGE

All models are supplied with 95 x 105 mm stage with integrated 16 x 50 mm X-Y mechanical stage



● JR.1051

JuniorLED

185 (h) x 140 (w) x 320 mm (d) | 2.3 kg

CONDENSER

Fixed condenser N.A. 0.65

ILLUMINATION

- All versions are supplied with adjustable transmitted and reflected 1W LED illumination
- All models have internal rechargeable batteries with an external 100-240 V battery 5V USB-C charger/mains adapter

PACKAGE CONTENT

- Supplied with USB-C type mains adapter-charger, three rechargeable AA batteries, an Allen key, dust cover, user manual
- All packed in a polystyrene box

MODELS	Head	Eyepiece 10x/16 mm	Achromatic objectives 4x/0.10, 10x/0.25 and S40x/0.65	Achromatic objectives 4x/0.10, 10x/0.25 and S60x/0.85	Achromatic objective S100x/1.25	illumination (incident & transmitted)	Rechargeable batteries
JR.1051	Monocular	•	•		o	2x 1 W LED	•
JR.1651	Monocular	•		•		2x 1 W LED	•

o = optional

ACCESSORIES

EYEPIECES

- JR.6010** Wide field eyepiece WF 10x/16 mm
- JR.6010-M** Wide field micrometric eyepiece WF 10x/16 mm
- JR.6016** Wide field eyepiece WF 16x/11 mm
- JR.6020** Wide field eyepiece WF 20x/11 mm
- JR.6025** Wide field eyepiece WF 25x/11 mm

OBJECTIVES

DIN Achromatic


- JR.7004** 4x/0.10 objective
- JR.7010** 10x/0.25 objective
- JR.7020** 20x/0.40 objective
- JR.7040** S40x/0.65 objective
- JR.7060** S60x/0.85 objective
- JR.7000** S100x/1.25 objective

MISCELLANEOUS

- JR.9992** Spare LED for JR.1051/1651
- JR.9992-E** Spare Electronics PCB for JR.1051/1651

DISPOSABLES

- PB.5155** Microscope slides 76 x 26 mm, ground edges, 50 pieces
- PB.5157-W** Microscope slides 76 x 26 mm, ground edges. White frosted side, 50 pieces per pack
- PB.5157-B** Microscope slides 76 x 26 mm, ground edges. Blue frosted side, 50 pieces per pack



- PB.5165** Cover glasses 18 x 18 mm, thickness 0.13-0.17 mm, 100 pieces
- PB.5168** Cover glasses 22 x 22 mm, thickness 0.13-0.17 mm, 100 pieces
- PB.5245** Lens cleaning paper, 100 sheets per pack
- PB.5255** Immersion oil, n = 1.515 (25 ml)
- PB.5274** Isopropyl alcohol 99% (200 ml)
- PB.5275** Cleaning kit: lens cleaning fluid, lint free lens tissue, brush, air blower, cotton swabs
- PB.5276** Microscope maintenance and servicing kit, 16 pcs: cleaning brush, 6 pcs screwdriver set, air blower, 3 pcs Allen key, 1.5, 2, 2.5 mm, lens cleaning fluid 20 ml, cleaning cloth 140 x 140 mm, 100 pcs lens tissue sheets, tube of maintenance grease, 10 ml bottle of oil, packed in a nice toolbox



● PB.5275



www.euromex.com