



S10/20 Series

S-20-LO





CLASSIC LINE - S10/20 Series

S-20-LO

Optical System	Greenough
Observation Tube	Binocular head
Inclination	45° inclined
Interpupillary distance	51-79mm
Diopter adjustment	On left tube, +/- 5 diopter
Eyepieces	Widefield WF10X/20mm
Objectives magnification	2X
Working distance	47mm
Stand type	Pole type stand
Head holder	For Ø14mm column with incident illumination
Focus mechanism	Coarse focusing system with tension adjustment
Focusing stroke	42mm
Incident illumination	Tungsten 12V/10W
Transformer	Internal
Power supply	220-240V (CE)
Accessories included	Dust cover, black/white stage plate
Dimensions LxWxH	170x115x240mm
Base	170x115mm
Column/Arm	Height 130mm
Weight	2,5Kg